

BWT AQU@SENSE MB

The Key to a Successful Contamination Control Strategy in Pharmaceutical Water

bwt-pharma.com

DBWT

AQU@Sense MB

TIME TO TAKE NEW PATHS

The AQU@Sense MB allows you to meet the requirements for continuous germ load monitoring in pharmaceutical water fully automatically.

"Out of the Box" Measurements

- All the solutions required for up to 1000 measurements are placed in a hermetically sealed cartridge,
- which takes just a few minutes to replace.
- Up to 99% of the cartridge is recycled.

Flow Cytometry

Flow cytometry is not just a high-precision method for counting all living microorganisms; it can do much, much more.

- It detects individual microorganisms, microorganisms in larger conglomerates, or bio film fragments and tells you the exact number of living cells.
- It creates a fingerprint of a microbiological population and can show any changes to it.

The Benefits

- Time saved for approval and troubleshooting
- Greater process reliability through trend analyses
- No "false positives" or "negative" results



Specific

The staining method provides maximum specificity for living organisms. No disruption caused by particles.

Online or Offline Fixed installation

or manual sampling

Fully Automatic

Up to 1000 measurements without any manual intervention



DBWT

Electronic Documentation

21 CFR Part 11-compliant electronic documentation

Proven Technology

AQU@Sense MB

Ø

The standard for drinking water and food for many years now

Primary Validation

In accordance with Ph.Eur. and USP requirements

FLOW CYTOMETRY

Quick, exact, reproducible and specific for intact cells







Precise identification of intact microorganisms by evaluating scattered light (blue line) and fluorescence signals (green and red lines)



The "DOT plot" compares the two fluorescence signals. The area under the gate shows a signal for each intact organism. This DOT plot is a kind of fingerprint for a population.

The bottom line: flow cytometry counts every single living cell



CONTAMINATION CONTROL STRATEGY

Annex 1 sets it out plainly: every system operator needs to have control over all their processes at all times, and must be able to prove that they are functioning correctly. That means seamless, riskbased monitoring has to be ensured.

The Three Levels of a CCS

The basis for all the other steps on the path towards establishing a comprehensive CCS is accurate knowledge of the actual situation in the system.

In a data-driven approach, that means:

- Recording knowledge of the current situation for all process steps promptly
- Carrying out an assessment based on a risk analysis
- Taking suitable measures (CAPA) to correct or prevent errors



hapter 5.1.6 and USP

The Plate Count Is Not the Solution

The conventional plate count (HPC) has many disadvantages, particularly regarding the time required to obtain measurement results. The use of an alternative method is recommended in Ph. Eur., Annex 1, provided that the information can deliver a scientifically sound measurement for assessing microbiological quality. The AQU@Sense MB meets these requirements.

Alternatives Are Already Available Today

A variety of technologies for determining biological quality quickly and accurately are already available on the market.

The AQU@Sense MB uses flow cytometry technology, which is probably the most advanced technology currently available and can also be used for online bioburden monitoring. It does not just show the presence of microbiological contamination, but also provides an exact cell count. Furthermore, the measurement includes additional information about the population of the microorganisms.



A typical measurement result of flow cytometry.

TECHNOLOGY YOU CAN TRUST

Primary validation has shown that flow cytometry significantly outperforms the conventional plate count in some areas under comparable conditions.

Flow cytometry is an established microbial counting technique. It has been used in laboratories and to assess drinking water for many years now. The AQU@Sense MB is the first flow cytometer that enables online monitoring of water quality in pharmaceutical systems. The basis for the validation of this application is a primary validation in accordance with the applicable regulations (Ph.Eur. 9.2, chapter 5.1.6 and USP 41, chapter 1223). This was carried out by BWT in collaboration with bNovate, Eindhoven University of Technology (TU/e) and the University of Applied Sciences Northwestern Switzerland (FHNW).

Range



not interior



Your Next Steps...

Validation for intended use

n the test r

Risk-benefit analysis

All species detected

Results

The AQU@Sense MB passed the primary validation in full!

HPC not infe

Parameter	Results
Accuracy	✓
Precision	✓
Specificity	✓
Quantification limit	✓
Linearity	✓
Measurement range	✓
Robustness	✓
Ruggedness	✓



ONE INSTRUMENT, MULTIPLE APPLICATIONS

The AQU@Sense MB can be used online or offline as a standalone device.

IN PROCESS CONTROL

Requirements for the method

- Specific
- Reproducible
- Quick
- Cost-effective

Benefits of AQU@Sense MB

- Keep a global overview and control each step of production
- Monitoring of all process steps for biological contamination, from the incoming water to the end product
- Early detection (trending) of signs of contamination with microorganisms to take preventative measures
- Reduced costs for water samples and associated tests
- Time saved during investigations

Trend Monitoring

The AQU@Sense MB lets you monitor microbiological contamination changes promptly, long before critical values are reached.

Sanitization Monitoring

- Sanitization effect immediately demonstrable
- Faster approval following maintenance



Trend Monitoring



Detecting changes in the system continuously and rapidly.





Sanitization Monitoring

Sanitization effect immediately demonstrable



QUALITY CONTROL UNIT

Requirements for the method

- Validated
- Specific
- Reproducible
- Correct

Benefits of AQU@Sense MB

- Immediate proof of quality
- Immediate intervention when problems occur
- Rapid analysis

Assess different results correctly

- Total bioburden background signal
- Detection and trend analysis of individual cells from conglomerates or bio film fragments.

SO TIME IS ON YOUR SIDE

Do you want to investigate yesterday's problem or shape tomorrow's success?

The AQU@Sense MB gives you the data you need, quickly, reliably and conclusively. The AQU@ Sense MB is not just a high-precision technology for counting all living microorganisms; it can do much, much more.

- It detects both individual organisms and organisms in larger conglomerates or biofilm fragments and tells you the exact number of living cells contained within them.
- It creates a fingerprint of a microbiological population in your system and can show any changes to it.

Automatic Digital Documentation

- No manual processes requiring time-consuming documentation.
- Analogue and digital interfaces to control systems
- Simple visualisation of results for further evaluations

Minimal Maintenance

- Menu-guided cartridge replacement in minutes following 1000 measurements
- No calibration





SERVICE AND SUPPORT

Minimal Maintenance

The AQU@Sense MB is designed in such a way that scheduled maintenance is required just twice per year. It is performed by qualified BWT staff as part of the fully documented AQU@Service.

Benefits

- No calibration required
- SST based on standard maintenance
- High availability as maintenance can be planned
- Transparency with maintenance and operating costs known in advance
- Hotline

Service Range

- Integration planning
- System integration
- Documentation
- Control system programming
- Automation adjustment
- Testing
- Project execution
- Change control
- Expertise/suggestions for improvement





DBWT

TECHNICAL DATA

General

Measuring principle	Flow cytometry
Degree of protection	IP65
Ambient temperature	+5 °C to +35 °C
Humidity	10 to 90% RH
Dimensions (W × D× H)	350 × 240 × 373 mm
Weight	14 kg

Connections

Outputs	4 digital outputs and 2 × 4 to 20 mA analogue outputs, freely configurable
Inputs	4 digital inputs
Sampler input	¼″ Swagelok tube fitting, male
Sampler output	¼″ Swagelok tube fitting, male
Energy supply	DC 18 V, 1.4 A, ax.20 W
Memory card	32 GB
Data export	USB and Ethernet

Measurements

AQU@Sense MB

.

90 µl
Continuous or manual
30 min. to 6 hrs
20 min.
Max. 1000 measurements
100 to 200 ml per min.
40 °C
Ozone up to 100 ppb Hot water up to 85 °C (without rupping measurement)



BWT AQUA AG

FOR YOU AND PLANET BLUE.